


Product Datasheet

Cytochrome c oxidase subunit 4 isoform 1, mitochondrial antibody (Biotin) (orb239877)

Description	Cytochrome c oxidase subunit 4 isoform 1, mitochondrial antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human Cytochrome c oxidase subunit 4 isoform 1, mitochondrial protein (23-169AA)
Target	COX4I1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P13073
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.



Western blot analysis of rat brain tissue...