


Product Datasheet

D-3-phosphoglycerate dehydrogenase antibody (Biotin) (orb239895)

Description	D-3-phosphoglycerate dehydrogenase antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human D-3-phosphoglycerate dehydrogenase protein (2-251AA)
Target	PHGDH
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	043175
Expiration Date	12 months from date of receipt.



Western blot analysis of (Lane 1): EC109...