



Product Datasheet

Endoplasmin antibody (orb239904)


Description	Endoplasmin antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, WB
Immunogen	Recombinant Human Endoplasmin protein (27-799AA)
Target	HSP90B1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P14625
Dilution Range	WB:1:1000-1:5000, IHC-P:1:20-1:200, IF:1:50-1:200
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human pl...



Immunofluorescence analysis of HeLa cell...



Western blot analysis of MCF-7 whole cel...