

Product Datasheet

Fructose-bisphosphate aldolase A antibody (orb239924)

Description Fructose-bisphosphate aldolase A antibody

Species/Host Rabbit

Reactivity Human, Rat

Conjugation Unconjugated

Tested Applications ELISA, IHC, WB

Immunogen Recombinant Human Fructose-bisphosphate aldolase A protein (2-364AA)

Target Hemoglobin

Preservatives Preservative: 0.03% Proclin 300.
Constituents: 50% Glycerol

Form/Appearance Liquid

Storage Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.

Note For research use only


Isotype IgG

Clonality Polyclonal

Uniprot ID **P04075**

Dilution Range WB:1:500-5000, IHC-P:1:20-200, IF:1:50-200, IP:1:200-2000


Expiration Date 12 months from date of receipt.



Immunohistochemistry of paraffin-embedde...



Western blot aanalysis of EC109 whole ce...



Immunohistochemistry of paraffin-embedde...