

Product Datasheet

Fructose-bisphosphate aldolase Ag antibody (FITC) (orb239926)

Description	Fructose-bisphosphate aldolase Ag antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Immunogen	Recombinant Human Fructose-bisphosphate aldolase A protein (2-364AA)
Target	α 1-acidglycoprotein
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P04075
Dilution Range	WB:1:500-5000, IHC-P:1:20-200, IP:1:200-2000
Expiration Date	12 months from date of receipt.