



Product Datasheet

Glycogen phosphorylase, brain form antibody (orb239956)


Description	Glycogen phosphorylase, brain form antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Recombinant Human Glycogen phosphorylase, brain form protein (6-843AA)
Target	PYGB
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P11216
Dilution Range	WB:1:500-1:5000, IHC-P:1:20-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of mouse skeletal ...



brain form antibody



Immunohistochemical staining of human pr...