



## Product Datasheet

# Nicotinate phosphoribosyltransferase antibody (orb240029)

|                            |                                                                          |
|----------------------------|--------------------------------------------------------------------------|
| <b>Description</b>         | Nicotinate phosphoribosyltransferase antibody                            |
| <b>Species/Host</b>        | Rabbit                                                                   |
| <b>Reactivity</b>          | Human, Mouse                                                             |
| <b>Conjugation</b>         | Unconjugated                                                             |
| <b>Tested Applications</b> | ELISA, WB                                                                |
| <b>Immunogen</b>           | Recombinant Human Nicotinate phosphoribosyltransferase protein (1-310AA) |
| <b>Target</b>              | NAPRT                                                                    |
| <b>Preservatives</b>       | Preservative: 0.1mol/L NaCl, 15mmol/L NaN <sub>3</sub> , PH7.2           |
| <b>Form/Appearance</b>     | Liquid                                                                   |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.            |
| <b>Note</b>                | For research use only                                                    |
| <b>Isotype</b>             | IgG                                                                      |
| <b>Clonality</b>           | Polyclonal                                                               |
| <b>Uniprot ID</b>          | <b>Q6XQN6</b>                                                            |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                          |



Western blot analysis of HepG2 whole cell lysate.



Western blot analysis of HepG2 whole cell lysate.