



Product Datasheet

O-phosphoseryl-tRNA(Sec) selenium transferase antibody (orb240037)


Description	O-phosphoseryl-tRNA(Sec) selenium transferase antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, WB
Immunogen	Recombinant Human O-phosphoseryl-tRNA(Sec) selenium transferase protein (5-451AA)
Target	SEPSECS
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9HD40
Dilution Range	WB:1:500-1:2000, IHC-P:1:20-1:200, IF:1:50-1:200
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of PC3 cells...



Immunohistochemical staining of human ki...



Western blot analysis of 293T whole cell...