

Product Datasheet

Phosphatidylserine decarboxylase proenzyme antibody (orb240077)

Description	Phosphatidylserine decarboxylase proenzyme antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant Human Phosphatidylserine decarboxylase proenzyme protein (1-375AA)
Target	PISD
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9UG56
Expiration Date	12 months from date of receipt.