

Product Datasheet

Steroid 17-alpha-hydroxylase/17, 20 lyase antibody (HRP) (orb240151)

Description	Steroid 17-alpha-hydroxylase/17, 20 lyase antibody (HRP)
Species/Host	Rabbit
Reactivity	Human
Conjugation	HRP
Tested Applications	ELISA
Immunogen	Recombinant Human Steroid 17-alpha-hydroxylase/17,20 lyase protein (24-508AA)
Target	CYP17A1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P05093
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.