

Product Datasheet

Sulfite oxidase, mitochondrial antibody (orb240152)

Description Sulfite oxidase, mitochondrial antibody

Species/Host Rabbit

Reactivity Human, Rat

Conjugation Unconjugated

Tested Applications ELISA, IHC, WB

Immunogen Recombinant Human Sulfite oxidase, mitochondrial protein (80-545AA)

Target SUOX

Preservatives Preservative: 0.03% Proclin 300.
Constituents: 50% Glycerol

Form/Appearance Liquid

Storage Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.

Note For research use only


Isotype IgG

Clonality Polyclonal


Uniprot ID **P51687**

Dilution Range WB:1:1000-1:5000, IHC-P:1:20-1:200


Expiration Date 12 months from date of receipt.



Western blot analysis
of Hela whole cell...



mitochondrial
antibody



Immunohistochemical
staining of human
ki...