

Product Datasheet

Biogenesis of lysosome-related organelles complex 1 subunit 1 antibody (FITC) (orb240225)

Description	Biogenesis of lysosome-related organelles complex 1 subunit 1 antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Immunogen	Recombinant Human Biogenesis of lysosome-related organelles complex 1 subunit 1 protein (1-153AA)
Target	BLOC1S1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P78537
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.