



Product Datasheet

Galectin-9 antibody (orb240262)


Description	Galectin-9 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, WB
Immunogen	Recombinant Human Galectin-9 protein (1-323AA)
Target	LGALS9
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	O00182
Dilution Range	WB:1:500-1:5000, IHC-P:1:200-1:500, IF:1:50-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of mouse heart tis...



Immunohistochemical staining of human ad...



Immunofluorescence analysis of NIH/3T3 C...