

Product Datasheet

Tyrosine-protein kinase transR ROR1 antibody (Biotin) (orb240361)

Description	Tyrosine-protein kinase transR ROR1 antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human Inactive tyrosine-protein kinase transmembrane receptor ROR1 protein (30-391AA)
Target	ROR1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q01973
Expiration Date	12 months from date of receipt.