



Product Datasheet

Ubiquitin-conjugating enzyme E2 K antibody (orb240368)

Description	Ubiquitin-conjugating enzyme E2 K antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC
Immunogen	Recombinant Human Ubiquitin-conjugating enzyme E2 K protein (2-200AA)
Target	UBE2K
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P61086
Dilution Range	IHC-P:1:10-1:100, IF:1:50-1:200
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human pl...



Immunofluorescence analysis of HeLa cell...