

## Product Datasheet

### Arabidopsis thaliana L-ascorbate peroxidase 2, cytosolic antibody (FITC) (orb240392)

|                        |                                                                                        |
|------------------------|----------------------------------------------------------------------------------------|
| <b>Description</b>     | Arabidopsis thaliana L-ascorbate peroxidase 2, cytosolic antibody (FITC)               |
| <b>Species/Host</b>    | Rabbit                                                                                 |
| <b>Reactivity</b>      | A. thaliana                                                                            |
| <b>Conjugation</b>     | FITC                                                                                   |
| <b>Immunogen</b>       | Recombinant Arabidopsis thaliana L-ascorbate peroxidase 2, cytosolic protein (4-250AA) |
| <b>Target</b>          | APX2                                                                                   |
| <b>Preservatives</b>   | Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol                            |
| <b>Form/Appearance</b> | Liquid                                                                                 |
| <b>Storage</b>         | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.                          |
| <b>Note</b>            | For research use only                                                                  |
| <b>Isotype</b>         | IgG                                                                                    |
| <b>Clonality</b>       | Polyclonal                                                                             |
| <b>Uniprot ID</b>      | <b>Q1PER6</b>                                                                          |
| <b>Expiration Date</b> | 12 months from date of receipt.                                                        |