

Product Datasheet

Arabidopsis thaliana Protein ABSCISIC ACID- INSENSITIVE 5 antibody (FITC) (orb240404)

Description	Arabidopsis thaliana Protein ABSCISIC ACID-INSENSITIVE 5 antibody (FITC)
Species/Host	Rabbit
Reactivity	A. thaliana
Conjugation	FITC
Immunogen	Recombinant Arabidopsis thaliana Protein ABSCISIC ACID- INSENSITIVE 5 protein (1- 442AA)
Target	ABI5
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or - 80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9SJN0
Expiration Date	12 months from date of receipt.