

Product Datasheet

Arabidopsis thaliana Trihelix transcription factor GT-1 antibody (HRP) (orb240409)

Description	Arabidopsis thaliana Trihelix transcription factor GT-1 antibody (HRP)
Species/Host	Rabbit
Reactivity	A. thaliana
Conjugation	HRP
Tested Applications	ELISA
Immunogen	Recombinant Arabidopsis thaliana Trihelix transcription factor GT-1 protein (1-406AA)
Target	GT-1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9FX53
Expiration Date	12 months from date of receipt.