

Product Datasheet

subtilis Alpha-galactosidase antibody (Biotin) (orb240419)

Description	subtilis Alpha-galactosidase antibody (Biotin)
Species/Host	Rabbit
Reactivity	Bacteria
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Bacillus subtilis Alpha-galactosidase protein (1-432AA)
Target	mela
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	O34645
Expiration Date	12 months from date of receipt.