

Product Datasheet

subtilis Inositol 2-dehydrogenase/D-chiro-inositol 3-dehydrogenase antibody (HRP) (orb240425)

Description	subtilis Inositol 2-dehydrogenase/D-chiro-inositol 3-dehydrogenase antibody (HRP)
Species/Host	Rabbit
Reactivity	Bacteria
Conjugation	HRP
Tested Applications	ELISA
Immunogen	Recombinant Bacillus subtilis Inositol 2-dehydrogenase/D-chiro-inositol 3-dehydrogenase protein (1-344AA)
Target	iolG
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P26935
Expiration Date	12 months from date of receipt.