

Product Datasheet

Baker's Pyruvate carboxylase 1 antibody (Biotin) (orb240451)

Description	Baker's Pyruvate carboxylase 1 antibody (Biotin)
Species/Host	Rabbit
Reactivity	Yeast
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Baker's yeast Pyruvate carboxylase 1 protein (18-470AA)
Target	PYC1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P11154
Expiration Date	12 months from date of receipt.