

Product Datasheet

Bordetella pertussis Serotype 2 fimbriasubunit antibody (orb240458)

Description	Bordetella pertussis Serotype 2 fimbriasubunit antibody
Species/Host	Rabbit
Reactivity	Bacteria
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant Bordetella pertussis Serotype 2 fimbrial subunit protein (27-207AA)
Target	fim2
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P05788
Expiration Date	12 months from date of receipt.