


Product Datasheet

Adiponectin antibody (orb240471)

Description	Adiponectin antibody
Species/Host	Rabbit
Reactivity	Bovine
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Bovine Adiponectin protein (18-240AA)
Target	ADIPOQ
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q3Y5Z3
Expiration Date	12 months from date of receipt.



Western blot analysis of recombinant Adiponectin protein (18-240AA)