

Product Datasheet

Vesicle-associated protein 7 antibody (Biotin) (orb240515)

Description	Vesicle-associated protein 7 antibody (Biotin)
Species/Host	Rabbit
Reactivity	Gallus
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Chicken Vesicle-associated membrane protein 7 protein (2-181AA)
Target	VAMP7
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q5ZL74
Expiration Date	12 months from date of receipt.