

Product Datasheet

Ornithine carbamoyltransferase chain F antibody (Biotin) (orb240535)

Description	Ornithine carbamoyltransferase chain F antibody (Biotin)
Species/Host	Rabbit
Reactivity	E. coli
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Escherichia coli Ornithine carbamoyltransferase subunit F protein (2-334AA)
Target	argF
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P06960
Expiration Date	12 months from date of receipt.