

Product Datasheet

Aerobic glycerol-3-phosphate dehydrogenase antibody (Biotin) (orb240557)

Description	Aerobic glycerol-3-phosphate dehydrogenase antibody (Biotin)
Species/Host	Rabbit
Reactivity	E. coli
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Escherichia coli Aerobic glycerol-3-phosphate dehydrogenase protein (1-501AA)
Target	glpD
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P13035
Expiration Date	12 months from date of receipt.