


Product Datasheet

ATP-dependent Clp protease adapter protein ClpS antibody (orb240560)

Description	ATP-dependent Clp protease adapter protein ClpS antibody
Species/Host	Rabbit
Reactivity	E. coli
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Escherichia coli ATP-dependent Clp protease adapter protein ClpS protein (1-105AA)
Target	clpS
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P0A8Q6
Expiration Date	12 months from date of receipt.



Western blot analysis of DH5 whole cell ...