


Product Datasheet

Ethanolamine ammonia-lyase light chain antibody (orb240576)

Description	Ethanolamine ammonia-lyase light chain antibody
Species/Host	Rabbit
Reactivity	E. coli
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Escherichia coli Ethanolamine ammonia-lyase light chain protein (1-295AA)
Target	eutC
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P19636
Expiration Date	12 months from date of receipt.



Western blot analysis of Rosetta bacteria...