

Product Datasheet

Mannose-6-phosphate isomerase antibody (FITC) (orb240606)

Description	Mannose-6-phosphate isomerase antibody (FITC)
Species/Host	Rabbit
Reactivity	E. coli
Conjugation	FITC
Immunogen	Recombinant Escherichia coli Mannose-6-phosphate isomerase protein (1-391AA)
Target	manA
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P00946
Expiration Date	12 months from date of receipt.