

Product Datasheet

Purine nucleoside phosphorylase DeoD-type antibody (Biotin) (orb240618)

Description	Purine nucleoside phosphorylase DeoD-type antibody (Biotin)
Species/Host	Rabbit
Reactivity	E. coli
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Escherichia coli Purine nucleoside phosphorylase DeoD-type protein (2-239AA)
Target	deoD
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or - 80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P0ABP8
Expiration Date	12 months from date of receipt.