



Product Datasheet

Ribokinase antibody (orb240621)

Description	Ribokinase antibody
Species/Host	Rabbit
Reactivity	E. coli
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Escherichia coli Ribokinase protein (1-309AA)
Target	rbsK
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P0A9J6
Expiration Date	12 months from date of receipt.



Western blot
analysis of
Rosetta
bacter...



Western blot
analysis of
DH5a whole
cel...