

Product Datasheet

Ribonucleotide reductase, class III, anaerobic antibody (FITC) (orb240627)

Description	Ribonucleotide reductase, class III, anaerobic antibody (FITC)
Species/Host	Rabbit
Reactivity	E. coli
Conjugation	FITC
Immunogen	Recombinant Escherichia coli Anaerobic ribonucleoside-triphosphate reductase protein (1-712AA)
Target	nrdD
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P28903
Expiration Date	12 months from date of receipt.