

Product Datasheet

Interleukin-8 antibody (HRP) (orb240662)

Description	Interleukin-8 antibody (HRP)
Species/Host	Rabbit
Reactivity	Gallus
Conjugation	HRP
Tested Applications	ELISA
Immunogen	Recombinant Gallus gallus Interleukin-8 protein (17-102AA)
Target	CXCL8
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P08317
Expiration Date	12 months from date of receipt.