


Product Datasheet

c-Jun complex subunit sigma-3 antibody (Biotin) (orb240700)

Description	c-Jun complex subunit sigma-3 antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA, IHC
Immunogen	Recombinant Human AP-1 complex subunit sigma-3 protein (1-104AA)
Target	Arginine esterase
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q96PC3
Dilution Range	IHC-P:1:20-200
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human te...