

Product Datasheet

Heat shock protein 75 kDa, mitochondrial antibody (orb240736)

Description	Heat shock protein 75 kDa, mitochondrial antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IP, WB
Immunogen	Recombinant Human Heat shock protein 75 kDa, mitochondrial protein (60-308AA)
Target	TRAP1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q12931
Dilution Range	WB:1:1000-1:5000, IF:1:50-1:200, IP:1:200-1:2000
Expiration Date	12 months from date of receipt.