

Product Datasheet

Protein phosphatase 1 regulatory subunit 11 antibody (FITC) (orb240746)

Description	Protein phosphatase 1 regulatory subunit 11 antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Immunogen	Recombinant Human E3 ubiquitin-protein ligase PPP1R11 protein (1-126AA)
Target	PPP1R11
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	O60927
Expiration Date	12 months from date of receipt.