

Product Datasheet

Klebsiella pneumoniae Beta-lactamase NDM-1 antibody (Biotin) (orb240775)

Description	Klebsiella pneumoniae Beta-lactamase NDM-1 antibody (Biotin)
Species/Host	Rabbit
Reactivity	Bacteria
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Klebsiella pneumoniae Metallo-beta-lactamase type 2 protein (29-270AA)
Target	blaNDM-1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	C7C422
Expiration Date	12 months from date of receipt.