

Product Datasheet

Leukocyte immunoglobulin-like R subfamily B member 3 antibody (orb240786)

Description	Leukocyte immunoglobulin-like R subfamily B member 3 antibody
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant Mouse Leukocyte immunoglobulin-like receptor subfamily B member 3 protein (664-841AA)
Target	Lilrb3
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P97484
Expiration Date	12 months from date of receipt.