


Product Datasheet

Acetylcholine R subunit alpha antibody (Biotin) (orb240803)

Description	Acetylcholine R subunit alpha antibody (Biotin)
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Mouse Acetylcholine receptor subunit alpha protein (21-230AA)
Target	Chrna1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P04756
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.



Western blot analysis of HeLa whole lysa...