


Product Datasheet

Apolipoprotein E antibody (orb240838)

Description	Apolipoprotein E antibody
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Mouse Apolipoprotein E protein (19-311AA)
Target	Probable terminase, ATPase subunit
Preservatives	Preservative: 0.1mol/L NaCl, 15mmol/L NaN ₃ , PH7.2
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P08226
Dilution Range	WB:1:1000-5000
Expiration Date	12 months from date of receipt.



Western blot analysis of mouse lung(Lane...