

Product Datasheet

Beta-defensin 14 antibody (Biotin) (orb240855)

Description	Beta-defensin 14 antibody (Biotin)
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Mouse Beta-defensin 14 protein (23-67AA)
Target	Defb14
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q7TNV9
Expiration Date	12 months from date of receipt.