

Product Datasheet

C-C motif chemokine 9 antibody (orb240878)

Description	C-C motif chemokine 9 antibody
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant Mouse C-C motif chemokine 9 protein (22-114AA)
Target	Ccl9
Preservatives	Preservative: 0.1mol/L NaCl, 15mmol/L NaN ₃ , PH7.2
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P51670
Expiration Date	12 months from date of receipt.