


Product Datasheet

CD82 antibody (Biotin) (orb240887)

Description	CD82 antibody (Biotin)
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Mouse CD82 antigen protein (111-227AA)
Target	Cd82
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P40237
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.



Western blot analysis of (Lane 1):jurkat...