


Product Datasheet

Mast Cell Chymase antibody (orb240890)

Description	Mast Cell Chymase antibody
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant Mouse Chymase protein (22-246AA)
Target	Cma1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P21844
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.



Western blot analysis of all lane lysate...