


Product Datasheet

Egl nine homolog 3 antibody (orb240933)

Description	Egl nine homolog 3 antibody
Species/Host	Rabbit
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Mouse Egl nine homolog 3 protein (2-239AA)
Target	Egln3
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q91UZ4
Expiration Date	12 months from date of receipt.



Western blot analysis of Rat brain tissue...