



Product Datasheet

Guanine nucleotide-binding protein G(o) subunit alpha antibody (orb240953)

Description	Guanine nucleotide-binding protein G(o) subunit alpha antibody
Species/Host	Rabbit
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Mouse Guanine nucleotide-binding protein G(o) subunit alpha protein (2-354AA)
Target	Gnao1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P18872
Dilution Range	WB:1:500-2000
Expiration Date	12 months from date of receipt.



Western blot analysis of rat brain tissue...



Western blot analysis of Rat brain tissue...