

Product Datasheet

Hemoglobin subunit beta-2 antibody (HRP) (orb240971)

Description	Hemoglobin subunit beta-2 antibody (HRP)
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	HRP
Tested Applications	ELISA
Immunogen	Recombinant Mouse Hemoglobin subunit beta-2 protein (2-147AA)
Target	Hbb-b2
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P02089
Expiration Date	12 months from date of receipt.