


Product Datasheet

Hepatitis A virus cellular R 1 homolog antibody (orb240972)

Description	Hepatitis A virus cellular R 1 homolog antibody
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Immunogen	Recombinant Mouse Hepatitis A virus cellular receptor 1 homolog protein (22-237AA)
Target	Havcr1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q5QNS5
Dilution Range	IHC-P:1:20-1:200
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human ki...