

Product Datasheet

Hyaluronan synthase 2 antibody (orb240980)

Description	Hyaluronan synthase 2 antibody
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant Mouse Hyaluronan synthase 2 protein (67-374AA)
Target	Has2
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P70312
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.