

Product Datasheet

Interferon epsilon antibody (orb240996)

Description	Interferon epsilon antibody
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant Mouse Interferon epsilon protein (22-192AA)
Target	Ifne
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q80ZF2
Expiration Date	12 months from date of receipt.